pempa®

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CE 0197

EN

Chapter 1	Introduction	1
Chapter 2	Setting the Meter	
2.1	Installing/ Replacing the Batteries.....	2
2.2	Set the Date and Time	2
2.3	Code the Meter.....	2
Chapter 3	How to Perform a Test	
3.1	Perform The Test	3-5
3.2	Alternate Site Testing (AST)	5-6
Chapter 4	Meter Memory Function	7
Chapter 5	Control Solution Test	8
Chapter 6	Care and Maintenance	
6.1	Storing Your Meter and Strip	9
6.2	Cleaning and Caring for Your Meter	9
Chapter 7	Error Message and Trouble Shooting	10-11
Chapter 8	Specification	12

Chapter 1 Introduction

Introduction

Please read carefully before using meter kit, and consult healthcare professional before making any important medical decision. Please contact your local customer service for further assistance with the product.

Normal Environmental Conditions

This meter designed under the following conditions:

- Indoor use
- Overvoltage category II
- Pollution degree 2

Electromagnetic Compatibility

This meter meets the electromagnetic compatibility, emission, and immunity, and the requirements of IEC 61326-2-6, ISO 18113-5 and IEC 61010-1.

Intended Use:

- *In vitro* diagnostic use only.
- Measuring blood glucose in fresh capillary whole blood from fingertip, palm, or forearm.
- Measuring total cholesterol/ uric acid in fresh capillary whole blood from fingertip.
- The meter can be used by laypersons or healthcare professionals.

Principles of the Examination Method:

- Electrochemical biosensor technology.


The meter is plasma-calibrated by reference instruments, which are traceable to the following standard reference materials and methods.

Test	Standard	Method
Glucose	NIST SRM 917	Glucose Oxidase
Total Cholesterol	NIST SRM 911	Abell / Kendall
Uric Acid	NIST SRM 913	Uricase / UV

Chapter 2 Setting the Meter

2.1 Installing/ Replacing the Batteries

This meter uses a CR2032 battery. Please remove the plastic tab under the battery before using. **Note:** Dispose the batteries according to your local environmental regulations.

2.2 Set the Date and Time

Press "S" Button (3 secs) → One "Beep" Sound → Setting Mode → Turn Off Automatically after Setting

- Setting order: Year/ Month/ Date/ Hour/ Minute
- Press "M" button to advance one unit, "S" button to enter next setting.


Note: Correct setting is important while managing your health records.

2.3 Code the Meter

- Code your meter when you first use it or open a new vial of strips.
- Make sure the meter is off before you insert the code strip.

- Make sure the codes on screen, code strip, and strip vial label are the same.

Note: With Auto Strip Recognition function, once you code your meter, the meter will switch to the test mode automatically when you insert a strip.


Chapter 3 How to Perform a Test

Materials you need to perform a test:

PEMPA Meter/ Test Strip/ Lancing Device/ Lancets/ Tissue or Cotton Ball with 75% Ethanol or Disinfection Wipes

3.1 Perform the Test

- Wash and clean your hands with disinfection wipes, and make sure your hands are dry before testing.



Insert lancet firmly.

Remove and save the protective cap.

Recap lancing device and adjust penetration depth.

Take a strip and close the vial immediately.

Insert the strip.

- 6  Make sure the code number is correct.
- 7  Pull the sliding barrel till it clicks.
- 8  Press the release button to sample.
- 9  Wipe off first drop of blood and start sampling.
- 10  Touch the blood sample with strip.

Glucose/ Uric Acid/
Total Cholesterol




Fill up ok

Fill up FAILED

Note: Insufficient fill up of blood sample could lead to inaccurate or failed test result.
Do not refill the test strip.

- 11  Result will show after countdown. Then discard the used strip to biohazard container.
- 12  The meter will turn off automatically.
- 13  Pierce used lancing into protective cap.
- 14  Discard used lancing to biohazard container.
- 15  Recap lancing device and storage.

Note:

- Please finish the test within 5 minutes or the meter will turn off automatically.
- The meter will not turn on if you insert the wrong end or wrong side of strip.
- Marked open date on new open strip vial, do not use expired strip.
- Dropping, bumping or other violent impact will damage the meter or cause malfunction.
- Do not use the meter in an environment with possible magnetic, electromagnetic, and radioactive interferences.

**Warning:**

- Please follow local regulations to discard used test strips and lancets.
- Used test strips, lancets and any other material that has been in contact with blood should be treated as potential biohazards.
- If user has infectious disease, the used test materials could be sources of infection.
- Lancets cannot be reused.
- Always use certified lancets to ensure safety.

3.2 Alternate Site Testing (AST)

You can test your glucose from fingertip, palm or forearm. Taking blood from palm or forearm could reduce the pain, but the glucose level changes faster. These differences may cause wrong medical decision.

Note: Please consult healthcare professional before AST sampling.

Suitable timing to acquire blood sample from alternate sites:

- Routinely before meal.
- Prior or 2 hours after meal/ short-acting or rapid-acting insulin analogue/ exercise.

DO NOT test from alternate sites:

- During or less than 2 hours after meal/ short-acting or rapid-acting insulin analogue/ exercise.
- When you think your glucose level is low or unaware of your low blood glucose condition.
- When you are examined for hypoglycemia or hyperglycemia.
- Your AST test result does not match your health condition.
- When you are ill, or you are operating machinery or driving a car.

Palm sampling

- No visible veins.
- Away from deep palm prints.

**Forearm sampling**

- Away from bones, visible veins and hair.

**Sampling from an Alternative Site:**

1. Repeat the steps 1-7 in Chapter 3.1. (Replace the lancing device tip with adjustable AST tip.)
2. Hold the lancing device against sampling site, and press the release button.
Keep holding the lancing device against sampling site until sufficient sample formed.
3. Then repeat steps 10-15 in Chapter 3.1.

Note:

- Sampling from fingertip if your AST test result does not match your health condition.
- Repeat puncturing the same spot may cause soreness and calluses.
- Do not squeeze the site excessively. It may take longer for sufficient blood sample to form.
- Do not use smeared blood sample, please acquire new blood sample.
- If you continue failing in getting enough blood samples, please try to get lancets in lower gauge or sample from fingertip instead.

Chapter 4 Meter Memory Function

- Glucose - 360 results, capable in counting 7-, 14-, 21- and 28- days average
- Total Cholesterol - 50 results
- Uric Acid - 50 results

The latest test result will replace the oldest when the records exceed maximum memory capacity. The memories start record from M1 to M360 or M1 to M50.

Directions for Checking Memories:

No Test Strip in the Meter → Press "M" Button (3 secs) → A Short "Beep" → Enter the Memory Mode → Full Display on Screen → Press "M" Button to Switch Mode (GLUC/ CHOL/ UA) → Press "S" Button to Confirm → Press "S" Button for Next Test Record/ Press "M" Button for Previous Record → Press "M" Button (3 secs) to Turn Off

- In GLU memory mode, it will display 7-, 14-, 21-, 28- days average first.
- Once you enter one memory mode, you cannot switch. You need to turn off the meter (Press "M" button for 3 seconds) and enter the memory mode again.

Memory Records are shown as follows:

Chapter 5 Control Solution Test


Control solution is used to check the performance of the kit.

The meter kit should be check:

- When the meter and strip do not work properly.
- When the test result is unusual or inconsistent.

Control solution range is shown as follow:

Please check your strip vial label for exact range.

Perform a Control Test:

Insert an Unused Strip → Make Sure the Code Number is Correct → Shake the Control Solution Well → Discard First Three Drops → Put One Drop onto a Clean Surface → Touch the Control Solution with Strip Sample Inlet → Test Result Display after Count Down → Compare the Result with the Range Listed on Test Strip Vial

Note:

- Do not reuse the test strip.
- Marked open date on new open control solution.
- Do not use the meter if the control test is out of range.
- If the control test keeps result wrongly, please contact an authorized distributor.
- Control solutions are not included. Please contact an authorized distributor for purchasing.

Chapter 6 Care and Maintenance

6.1 Storing Your Meter and Strip

Meter:

- Avoid bump or violent behavior.
- Do not use in extremely dry environment. It may cause static discharges.
- Do not use under electromagnetic radiation, ex: electrical equipment.
- Do not disassemble the meter for any reason.
- Keep meter clean by wiping the exterior appearance with tissues or lint-free cloth.

Strip:

- Do not store in high humidity environment, or expose directly to sunlight.
- Do not freeze or refrigerate the meter and strips.
- Keep your hand dry and clean while handling the strips and performing the test.

6.2 Cleaning and Caring for Your Meter


Please use soft cloth slightly damp with one of the following solution to clean meter surface.

- 75% alcohol
- Super Sani-Cloth disposable wipes
- Mild dishwashing liquid with water
- 10% household bleach solution and 90% water

Note:

- Do not allow any other wet cloth or liquid.
- Do not allow any liquid run in or around the test port and battery cover.
- Make sure the meter is completely dry before use.
- Protection impairment if used in a manner not specified by the manufacturer.

Chapter 7 Error Message and Trouble Shooting

Message	Cause	Solution
	<ul style="list-style-type: none">• Problem with code strip.• Problem with test strip.• Insert strip improperly.	Repeat the coding procedure (Chapter 2.3) and insert the strip again. If the problem persists, please contact local distributor for service.
	<ul style="list-style-type: none">• Low battery.	Replace with new battery.
	<ul style="list-style-type: none">• Problem with code strip or meter.• Insert code strip improperly.	Repeat the coding procedure (Chapter 2.3). If the problem persists, please contact local distributor for service.
	<ul style="list-style-type: none">• Incorrect meter operating temperature.	Repeat the test after meter return to operating temperature. If the problem persists, please contact local distributor for service.
	<ul style="list-style-type: none">• Used strip.• Damped strip.	Follow Chapter 3.1 and repeat the test with a new strip. If the problem persists, please contact local distributor for service.

E-9	<ul style="list-style-type: none"> • Test incomplete due to removing the strip during measuring. 	Follow Chapter 3.1 and repeat the test with a new strip. Do not remove the strip before the test is completed.
E-A	<ul style="list-style-type: none"> • Problem with strip. 	Follow Chapter 3.1 and repeat the test with the same strip. If the problem persists, please contact local distributor for service.
E-1	<ul style="list-style-type: none"> • Improper code strip. 	Repeat the coding procedure (Chapter 2.3). If the problem persists, please contact local distributor for service.
E-8	<ul style="list-style-type: none"> • Sample volume not enough. 	Repeat the test with a new strip, and make sure the sample volume is enough. If the problem persists, please contact local distributor for service.
HI	<ul style="list-style-type: none"> • Test result is higher than the range listed on Chapter 8. 	Follow Chapter 3.1 and repeat the test with the same strip. If the problem persists, please contact local distributor for service.
Lo	<ul style="list-style-type: none"> • Test result is lower than the range listed on Chapter 8. 	Follow Chapter 3.1 and repeat the test with the same strip. If the problem persists, please contact local distributor for service.

Chapter 8 Specification

Test Sample	Fresh Capillary Whole Blood
Measuring Time	GLU: 5 seconds; CHOL: 26 seconds; UA: 15 seconds
Measuring Range	GLU: 20-600 mg/dL (1.1-33.3 mmol/L); CHOL: 100-400 mg/dL (2.59-10.35 mmol/L); UA: 3-20 mg/dL (0.18-1.19 mmol/L)
Sample Volume	GLU: 0.7 µL; CHOL: 0.8 µL; UA: 1 µL
Storage & Transportation Condition	4-30°C (39-86°F)
Operation Temperature	10-40°C (50-104°F)
Storage & Transportation Relative Humidity	10-90%
Memory	460 Test Results (GLU: 360 ; CHOL: 50 ; UA: 50)
Battery Type	One 3V (CR2032) lithium battery
Battery Life	Approximately 1,000 tests
Dimensions	86*57*17 mm
Weight	About 48g (with battery)
Altitude	10,000 feet (3048 m)
Expected Service Life	5 years

Note:

- Please refer to the strip insert for accuracy, precision, limitation, and other important information.

Contents of the Kit (please check the meter outer box for exact detail)

- PEMPA PLUS Meter (BK6-40M) (with a CR2032 3V battery)
- User's Manual
- Quick Guide
- Carry Bag
- Lancing Device

Optional (not included in the standard kit package, please contact your local distributor for ordering)

- PEMPA Glucose Test Strips (BK6-G)(with code card & Insert)
- PEMPA Total Cholesterol Test Strips (BK-C2)(with code card & Insert)
- PEMPA Uric Acid Test Strips (BK-U1)(with code card & Insert)
- Lancets (Please refer to package for manufacture information.)

Note: After purchasing, if the contents are damaged, please contact authorized distributor immediately.